

IMIXCube™


The IMIXCube, the first self supporting overhead tube crane, provides the flexibility and many positioning advantages of an overhead tube crane without the cost and hassle of unistrut.



**Overhead flexibility
without Unistrut!**

we™
are
dr

IMIXCube™



Two different versions

The Installation cube is available in two different versions:

Version 1 is intended for really small rooms requiring only 2,5m (8'3") ceiling height.

Version 2 gives you an increased stroke in the z-direction for rooms with ceiling height of 2,8m (9'3") or more.

Usable stroke

The length of the stroke of the ceiling suspended x-ray tube:

X: 1850mm (6'1")

Y: 2750mm (9'1")

Z1: 1450mm (4'9"), Version 1

Z2: 1700mm (5'7"). Version 2

Focuspoint to floor¹

approx. 300–1750mm (1'–5'9"), Version 1

approx. 300–2000mm (1'–6'7"), Version 2

¹) Depending on the selected tube.

* Specifications subject to change without notice



IMIX is committed to providing high-quality, value-driven DR solutions – everytime.

IMIX Americas Inc.

22560 Glenn Drive, Suite 116, Sterling, VA 20164
 PHONE +1 [703] 430 3111, PHONE 1-888-IMIX RAD
 FAX +1 [703] 430 8808, WEB imixadr.com

