

IMIXDetectors™

IMIXSystems


family of Digital Detectors combine high reliability with outstanding image quality. **IMIXSystems Detectors** are delivered complete with IMIX's comprehensive Control Station SW. The result is a ready-to-integrate platform providing unparalleled value in Digital Radiography.



IMIXDetectors™

IMIXinsight™ CCD Detector Specification

	FOCUS Series	PRO Series	ELITE Series
Number of Pixels	9 MP	16 MP	16 MP
Imaging area	16" x 16"	17" x 17"	17" x 17"
Image acquisition	16-bit	16-bit	16-bit
Spatial Resolution	> 3.1 lp/mm	> 4.2 lp/mm	> 4.2 lp/mm
Scintillator	Gd202S	Gd202S	CsI (ACS-High Light)
140 cm focus grid	60 line	80 line	80 line
3 field ion chamber	Optional	Optional	Optional



IMIX Control Station

- Image review within 10 sec of exposure
- Windows based image review tools
- 20" flat panel display
- Local storage for more than 1000 images
- DICOM 3.0 compliant
- Network connectivity

* Specifications subject to change without notice



IMIX is committed to providing high-quality, value-driven DR solutions – everytime.



IMIX Americas Inc.

22560 Glenn Drive, Suite 116, Sterling, VA 20164
 PHONE +1 [703] 430 3111, PHONE 1-888-IMIX RAD
 FAX +1 [703] 430 8808, WEB imixadr.com

