

IMIX Slimline™ Detectors

IMIX Slimline™ Detectors

provide high performance DR imaging in a slim package, perfect for dual detector, retrofit and other applications where space constraints are an issue.



IMIX Slimline™ Detectors

IMIX Slimline Detectors™ Specifications

	SLIMLINE	SLIMLINE T
Imaging Technology	a-Si / TFT	a-Si / TFT
Scintillator	CsI or DRZ Plus	CsI or DRZ Plus
Configuration	Fixed Bucky Mount	Tethered
Field of View	17" x 17"	14" x 17"
Spatial Resolution	3.6 lp/mm	3.6 lp/mm
Pixel Count	9 MP	7½ MP
Pixel Pitch	139 microns	139 microns
Pixel Matrix	3K x 3K	2½K x 3K
Image Acquisition	14 bit	14 bit
Warranty	1 year	1 year

Control Station

- Image review in 1 second
- Windows based image review tools
- 20" flat panel video display
- Local storage for more than 1'000 images
- DICOM 3.0 compliant
- Network connectivity



* Specifications subject to change without notice

IMIX Slimline T™

provides added flexibility to make existing x-ray systems cassette-FREE.



IMIX is committed to providing high-quality, value-driven DR solutions – everytime.

IMIX Americas Inc.

22560 Glenn Drive, Suite 116, Sterling, VA 20164
 PHONE +1 [703] 430 3111, PHONE 1-888-IMIX RAD
 FAX +1 [703] 430 8808, WEB imixadr.com

